HT-IMPLANT

User's Manual

"HT IMPLANT" Implant/Surgery System

HT-IMPLANT

Manufactured by	
	By Dental • A company of ACTEON Group Via Vecchia Provinciale Lucchese, 49 51030 SERRAVALLE PISTOIESE Italy
Distributed by	
	SATELEC • A company of ACTEON Group 17 av. Gustave Eiffel BP 30216 33708 MERIGNAC cedex FRANCE
	EU Directive 93/42/EEC has been used as a basis in the design and manufacture of this medical product. A product with this mark meets the requirements of the relevant EC directives, i.e. the applicable standards in Europe
IEC 60601-1 EN ISO 13485	This equipment is manufactured according to current regulations and standard (IEC 60601-1) and according to the EN ISO 13485 quality control certification system.

Please indicate the product serial number in all requests.

ENGLISH	4
FRANÇAIS	27
ESPAÑOL	52
DEUTSCH	77
ITALIANO	104
NEDERLANDS	130
PORTUGUÊSE	156
SVENSKA	182
ΕΛΛΗΝΙΚΑ	206

ENGLISH

CONTENTS

I - INTRODUCTION	6
II - WARNINGS.....	6
III - DESCRIPTION	9
3.1 PHYSICAL DESCRIPTION.....	9
3.2 TECHNICAL DESCRIPTION	9
a. LCD screen/control Keypad	9
b. Control unit back panel	9
c. Control footswitch.....	9
d. Technical characteristics.....	10
IV- INSTALLATION / FIRST USE	10
4.1 UNPACKAGING THE DEVICE	10
4.2 RECOMMENDATIONS.....	10
4.3 INSTALLATION	10
V- IRRIGATION SOLUTION	11
VI- IRRIGATION LINE	11
VII- ADJUSTMENTS / SETTINGS / INTERFACE.....	11
7.1 START UP	11
7.2 SETTINGS	11
7.3 TORQUE CALIBRATION.....	12
7.4 PROGRAM SETTINGS ADJUSTMENT.....	12
a. Contra-angle selection.....	12
b. Speed adjustment	12
c. Torque adjustment	12
d. Flow	13
e. Reverse.....	13
f. Motor.....	13
g. Pump	13

h. Footswitch	13
VIII- USING THE DEVICE.....	13
IX - SHUTTING DOWN THE DEVICE	14
X - ROUTINE MAINTENANCE / STERILIZATION	15
10.1 IRRIGATION LINE MAINTENANCE	15
10.2 CONTRA-ANGLE MAINTENANCE	15
10.3 ROTARY INSTRUMENT MAINTENANCE	15
10.4 DEVICE MAINTENANCE.....	15
10.5 MICRO-MOTOR AND CORD MAINTENANCE	16
10.6 MICRO-MOTOR HOLDER MAINTENANCE	16
XI - MONITORING / PREVENTIVE AND CORRECTIVE MAINTENANCE.....	16
11.1 MONITORING.....	16
11.2 PREVENTIVE AND CORRECTIVE MAINTENANCE.....	16
11.3 FUSE REPLACEMENT	17
11.4 OPERATION FAULTS	18
XII - ELECTROMAGNETIC COMPATIBILITY	21
12.1 ELECTROMAGNETIC EMISSIONS.....	21
12.2 ELECTROMAGNETIC IMMUNITY	22
12.3 RECOMMENDED SEPARATION DISTANCES	23
XIII - DISPOSAL AND RECYCLING.....	24
XIV - MANUFACTURER'S LIABILITY	24
XV - REGULATIONS	24
XVI - SYMBOLS	24
XVII - CUSTOMER RELATIONS	233

I - INTRODUCTION

Congratulations! You have just taken possession of your HT IMPLANT™, a highly-engineered drive unit for dental surgery applications! This powerful unit has a reliable electronic speed and torque control with an integrated coolant pump. All on-demand operations during the implantology session can be accomplished with complete safety and extreme accuracy.

- The brushless induction micro motor is able to reach a speed of 36,000 Revolutions Per Minute (RPM).
- The E-Type - ISO 3964-type coupling means that the motor can be used with most surgical contra-angles without internal light or spray on the market. Safe, powerful, reliable, silent, light weight and user friendly.
- The rotor of the pump is protected and the tube can be easily installed and/or changed. Our innovative "NO DROP" system design alleviates the problem of solution dropping through the handpiece once the micro motor has been halted. For the correct operation of the peristaltic pump, use only: tubings distributed by SATELEC.

In order to get the most out of this high-technology device and ensure it has a long working life, please read this manual carefully before using or maintaining it.

Sentences preceded by the symbol  are points to which you should pay particular attention.

II - WARNINGS

CAUTION:

United States Federal Law restricts the use of this device solely to qualified, trained and competent dental health practitioners or under their supervision.

To reduce the risks of accidents, the precautions stipulated below must be taken.

Device users:

- Use of the HT IMPLANT is restricted solely to qualified, trained and competent dental health practitioners in the normal context of their work.
- If you have received this device by error, please contact the supplier so that it can be removed.

Interactions:

- The device complies with applicable electromagnetic compatibility standards. The user should nevertheless ensure that any potential electromagnetic interference does not cause an additional risk (presence of radiofrequency emitters, electronic devices, etc.).

- Interference may occur when used on patients with cardiac pacemakers. This system emits electromagnetic fields, which means there are some potential risks. The malfunctioning of implantable devices such as cardiac pacemakers and ICDs (implantable cardioverter defibrillator) is possible:

- Ask patients and users if they have an implanted device before using this product. Explain the circumstances to them
- Weigh the risks and benefits and contact your patient's cardiologist or appropriate qualified healthcare professional prior to performing the treatment
- Keep this product away from implanted devices
- Make appropriate emergency provisions and take immediate action if patients become ill

- Symptom including a raised heartbeat, irregular pulse, and dizziness may signal problems with a cardiac pacemaker or ICD.

Electrical connection:

- Your device must be connected to the electric power supply by a certified dental installation technician.
- To avoid any risk of electric shock, this device must be connected to an electric supply with a protective ground.
- The electric supply to which the device is connected must comply with the standards in force in your country.

Using the device:

- Do not use the device if it appears to be damaged or faulty.
- Turn the device off before unplugging the power cord.
- To unplug the power cord, grip the cord plug and hold the wall socket.
- Never use any other irrigation solution containers than those intended for suspension from the supplied brackets.
- The device must only be used with bottles or bags of physiological saline.
- The capacity of the irrigation solution containers used must not exceed one liter.
- When the device is not to be used for a long period of time, unplug the device from the electric supply.
- Do not move the device during use.

	▪ Do not insert anything inside the motor ▪ The “Motor” is supplied Non-Sterile! ▪ Use only power cord provided (chapter III)
	▪ Use for intended purposes only. Failure to observe the operating instructions may result in the patient or user suffering serious injury. Before using this product, make sure that you have studied and understood the operating instructions. ▪ Do not install where there is a risk of an explosion. The Systems are not intended for operation in the presence of flammable anesthetics or gases. ▪ Do not disassemble or alter the System motor, console, or footswitch.
	▪ Connect mains cord to a properly grounded outlet only. ▪ Never touch drills, burs, or other handpiece tips when they are still rotating. ▪ Handpiece should only be attached when the motor has stopped running.

Environment:



- Do not cover the device or obstruct the ventilation vents.
- Do not immerse the device in liquid, and do not use it outdoors.
- Do not tilt the device at an angle greater than 5°.
- Do not place the device near a heat source.
- Make sure that the cords are not in a traffic path.
- The device should be stored in its original packaging in a safe place.
- The device is not designed for use in the presence of anesthetic gases or any other flammable gas.
- Do not expose the device to water vapor, or splashes.
- The device is not designed to work near ionizing radiation.
- Do not insert metal objects into the device (risk of electric shock, short-circuit or emission of hazardous substances).

Maintenance:

- Before and after each use, your device and accessories must be cleaned and disinfected with products recommended by the manufacturer.
- Do not use agents containing phenol, peracetic acid, peroxide and other agents that break up the oxygen.
- DO NOT (for any reason), lubricate the motor. After lubricating the handpiece, ensure that the lubricant does not penetrate into the motor.

Accessories:

- The device can only accept accessories designed by BY DENTAL or distributed by SATELEC for the particular use for which they are intended.
- Use of accessories from other manufacturers is a potential hazard for you or your patients.

Repair:

- Do not repair or modify the device without prior authorization from the manufacturer.
- In the case of a fault, contact the supplier of your device. Do not use unauthorized repairers, who might make your device dangerous for you and your patients.

If you have any doubt, contact an approved dealer or our customer support department:

- www.acteongroup.com
- E-mail: satelec@acteongroup.com

III - DESCRIPTION

3.1 PHYSICAL DESCRIPTION

- Electronic Control Console (Fig.1 - 1)
- Brushless Micro Motor (Fig.1 - 2)
- Footswitch (Fig.1 - 3)
- Motor Holder (Fig.1 - 4)
- Bracket (Fig.1 - 5)
- Peristaltic pump tube (Fig.1 - 6)
- Mains cord (Fig.1 - 7)

Note: All accessories are supplied non-sterile with the exception of tubing

3.2 TECHNICAL DESCRIPTION

a. LCD screen/control Keypad

The HT IMPLANT is adjusted by applying reasonable pressure to the keypad buttons (Fig.2).

- Motor ON/OFF (Fig.2-1)
- Pump ON/OFF (Fig.2-2)
- Pump Flow adjustment (Fig.2-3)
- Storage and quick access up to 5 different “Programs” (implant configurations) (Fig.2-3)
- Speed Increase/Decrease (Fig.2-4)
- Torque adjustment (Fig.2-5)
- Handpiece ratio selector (Fig.2-6)
- Forward/Reverse Selector (Fig.2-7)
- Calibrate handpiece (Fig.2-8)
- Motor Receptacle (Fig.2-9)
- Footswitch in (Fig.2-10)
- Bracket holder (Fig.2-11)

Information required by the user is displayed on the LCD screen.

b. Control unit back panel

1. Irrigation pump housing and irrigation tube insertion (Fig. 3- 1)
2. Mains cord in (Fig. 3- 2)
3. Main Power On/Off Switch (Fig. 3- 3)
4. Fuses and voltage selector (Fig. 3- 4)

c. Control footswitch

Motor variation to the pre-selected maximum speed (Fig. 4).

d. Technical characteristics

Manufacturer	By Dental - ACTEON Group - Pistoia - Italy	Nominal Power	70 VA
Distributor	SATELEC - ACTEON - France		
Model	HT IMPLANT	Fuses	2x1AT(230V) - 2x2AT(115V) 5 x 20
Dimensions (W x H x D)	240 x 100 x 210 mm	Protection	Class I
Materials	NORYL KB 250	Type	BF
Micro motor	36000 rpm / Torque: 28.6 mN.m	Class MD	Ila rule 9 in Annex IX 93/42/EEC
Noise	< 65 dBA	Console	IP21
Supply Voltage	230 VAC/115 VAC	Micro motor	IP65
Frequency	50 Hz /60 Hz	Footswitch	IP21
Not suitable for use in presence of flammable anesthetics or oxygen.			
Conditions of use: Temp +18°C/+40°C (+64°F/+104°F) RH < 80%. Altitude restriction: 3.000m (10.000 ft)			
Shipping and Storage conditions: Temp +5°C/+65°C (+41°F/+149°F) RH < 20-95% non condensing.			

IV- INSTALLATION / FIRST USE

4.1 UNPACKAGING THE DEVICE

Upon reception of the device, check for any damage caused in transit.
Contact your supplier if necessary.

4.2 RECOMMENDATIONS

Have your device connected to the power supply by a certified dental installation technician. The electrical connection of the HT IMPLANT must comply with the standards applicable in your country.
A grounded power supply socket must be used.

4.3 INSTALLATION

Important:

Do not place the HT IMPLANT close to or on top of another device.

Do not place the power cord and the footswitch cord in a wire cover or in a cable gland.

1. Remove the unit and its accessories from their original packaging and place on a flat surface. Not install near heat sources, direct or indirect.
2. Verify the contents of the box according to the section “3.1 PHYSICAL DESCRIPTION”.
3. Do not kink the motor hose since it may damage it.
4. All accessories are supplied non-sterile with the exception of tubing.
5. The serial number on the Unit's rear panel must be the same as the one on transport documents.
6. Allow at least 150 mm free space around the unit for cooling ventilation.

- Place the control unit in position on a fixed surface, horizontal or with a slope of no more than 5°.
- Check that the power switch is in position O (off).
- Connect the power cord to the receptacle of the device.
- Connect the power cord to a socket with a ground pin.
- Connect the control footswitch to the footswitch connector.
- Position the footswitch so that it is easily accessible for your feet.
- Insert the bracket into the bracket holder.
- Connect the micro-motor and cord to the connector.
- Place the micro-motor on its holder.
- Hang the bottle or bag of physiological saline from the bracket.

V- IRRIGATION SOLUTION

The HT IMPLANT is not designed to administer drugs. It must be used only with bottles or bags of physiological saline or sterile water with a capacity not exceeding one liter.

VI- IRRIGATION LINE

The HT IMPLANT must be used exclusively with disposable sterile irrigation lines distributed by SATELEC.

VII- ADJUSTMENTS / SETTINGS / INTERFACE

7.1 START UP

First switch on → unit screen (Fig. 5)

- The First row shows the Program in use [1 to 5], the torque in use and the motor sense.
- The Second row shows the irrigation flow rate [0-1-2-3], the ratio of the contra-angle in use and the speed in use.

In every program Programs [1 to 5] are set the same parameters by the manufacturer:

- Ratio: 20:1
- Speed: 1800 rpm
- Flow: 3 (Max)
- Torque: maxNcm
- Motor sense: FWD

The listed parameters are only examples. To prevent unnecessary risk, observe the guideline speeds given by the manufacturer of the rotating instruments and the guideline given by the manufacturer of the “implant System” in use.

7.2 SETTINGS

HT IMPLANT will store the operative parameters for each “program” as soon as they are set by the user. They will be maintained in the memory even if the central unit is switched off. For each “program” it is possible to set: speed, torque, ratio, direction of rotation and flow.

To reset the parameters of the manufacturer (Software unit reset), press the “PUMP ON/OFF” key while the unit is switched on (Fig. 6).

7.3 TORQUE CALIBRATION

Perform torque calibration, press “Calibration” key on front panel while holding the motor (not drilling). This function checks the handpiece mechanical inertia, in order to off-set this value. This operation will take a few seconds, the drill will run at different speed (Fig. 7).

- The drill must not be engaged in any operation when this function is performed.
- Calibration should be done before each surgical session or anytime the handpiece is changed or lubricated.

- *Do not perform calibration if the micro motor is not inserted in the receptacle on the front panel of the unit. In case of error, perform a “System Reset” as described previously.*

7.4 PROGRAM SETTINGS ADJUSTMENT

a. Contra-angle selection

Sets handpiece reduction ratio. See the available reduction ratio in the following table:

1:1	16:1	20:1	24:1	32:1	64:1	80:1
-----	------	------	------	------	------	------

b. Speed adjustment

SPEED +/-: Sets motor rotation speed (value in rpm x 20 steps). Speed values depend on contra-angle reduction rate. See the following table:

Reduction Rate	Max. RPM	Min. RPM
1:1	36,000	400
16:1	2250	25
20:1	1800	20
24:1	1500	16
32:1	1125	12
64:1	1125	6
80:1	450	5

Avoid excessive speed. Set the speed to a safe maximum for the application.

c. Torque adjustment

TORQUE +/-: increases or decreases the torque values (24 steps). Torque values are in Ncm and depend on contra-angle reduction rate. See the following table with the **handpiece 20:1**:

Ncm					
3.0	6.0	9.0	12.0	15.0	18.0
21.0	24.0	27.0	30.0	33.0	36.0
39.0	42.0	45.0	48.0	51.0	54.0
57.0	60.0	63.0	66.0	69.0	MAX

Torque accuracy is guarantee only with a contrangle correctly maintained. The torque values are only for handpieces that operate properly. Micro motor stops when set torque is reached. The torque value in counterclockwise rotation is always = Max Ncm

d. Flow

FLOW: Selects variable irrigation pump flow rate. See the following table:

0 = no irrigation	1 = MIN	2 = MED	3 = MAX
-------------------	---------	---------	---------

e. Reverse

REVERSE: Switch between forward and reverse motor rotation (when in reverse, console will emit a beeping tone).

f. Motor

MOTOR ON/OFF: Switches motor ON and OFF (if enabled).

g. Pump

PUMP ON/OFF: Switches pump ON and OFF (if enabled).

Always check the values before use. Danger of injury for incorrect use.

h. Footswitch

Motor control up to the pre-selected maximum speed.

VIII- USING THE DEVICE

Before turning on the central unit, check to see that you have completed the installation of ALL COMPONENTS.

Important:

- Do not disconnect the micro-motor cord when the device is switched on and the footswitch is pressed.
- Do not engage the instrument to or disengage it from the contra-angle when the micromotor is running.
- All accessories must be cleaned and disinfected before use.
- For your safety and that of your patient, the HT IMPLANT must not be used with accessories other than those supplied by BY DENTAL and/or SATELEC.
- Check the integrity of the device and its accessories before and after each use in order to detect any problem.
- If necessary, do not use the device and replace any faulty component.

The HT IMPLANT is prepared for use as follows:

- Make sure that the device is connected correctly and properly insulated.

- Micromotor connection:

- Connect the micro motor into the receptacle on the front of the central unit (Right - white - 8 Pins) (Fig. 2-10), (Fig.8)
- Connect the “E-type” ISO 3964 contra-angle to the micro motor and fit an instrument (drill, bur, etc.) to the contra-angle. (Fig. 9)
- Place the micro motor on the motor holder.

- Footswitch connection:

- Connect the footswitch into the receptacle on the front of the central unit (Left - Blue - 5 Pins) (Fig. 2-11)
- Place the footswitch on the ground.

While inserting the plug, align the guide on the connector in an upright position to fit into the corresponding female guide on the HT IMPLANT. Connectors have a spring retention system to help avoid accidentally disconnecting the unit's cords and cables.

- Irrigation connection:

- Hang the bottle or bag of physiological saline from the bracket.
- Insert the irrigation line to the pump (“IN” = perforator side ; “OUT” = contra-angle side) and insert the irrigation line perforating pin into the bottle or bag of physiological saline (Fig. 10).
- Fasten the irrigation line to the cord using line clips.
- Connect the end of the irrigation line to the contra-angle.

The manufacturer accepts no responsibility in the case of irrigation tubes other than those mentioned in this manual, being used.

- Switch on the device:

- Switch the power switch to position I (power on).
- Make any necessary setting adjustments (irrigation flow rate, speed, torque, program, etc.)
- Use the HT IMPLANT in accordance with good dental practices.

IX - SHUTTING DOWN THE DEVICE

Upon completion of the dental procedure:

- Remove the irrigation line perforating pin from the bottle or bag of physiological saline.
- Remove the bottle or bag from the bracket.
- Immerse the irrigation line perforating pin in a vessel containing distilled water.
- Rinse the irrigation line and the contra-angle by running the flush function until the recipient and irrigation lines are completely empty.
- Remove the irrigation line clips.
- Disconnect the irrigation line from the contra-angle and dispose of it in a safety container for soiled medical items.
- Remove the instrument fastened to the contra-angle.
- Disconnect the contra-angle from the HT IMPLANT micro-motor.
- Disconnect the micro-motor cord from the HT IMPLANT unit.
- Sw
itch the device off (O).

X - ROUTINE MAINTENANCE / STERILIZATION

Important:

The HT IMPLANT micro-motor and cord is not autoclavable.

Cleaning	Materials: Water, Neutral detergent, damp cloth, alcohol	To remove dirt or traces of drugs on the surface, wipe with a cloth slightly dampened with water, neutral detergent or alcohol where permitted. Be sure to remove any dirt and residue from all components.
Disinfection	Materials: Disinfectant tested and approved with demonstrable bactericidal, fungicidal and virucidal properties. Method: wipe with a cloth slightly dampened with disinfectant products.	
Drying	Materials: Soft cloth	Drying carefully with a soft and clean cloth all the components.
Maintenance	The micro motor requires no maintenance. - Check at least every 6 months the cable of the micro motor. In case of deterioration, it must be replaced at an authorized service center.	
Shipping and Storage	Temp +5°C/+65°C (+41°F/+149°F) RH < 20-95% non condensing.	

10.1 IRRIGATION LINE MAINTENANCE

The irrigation lines and clips distributed by SATELEC are disposable and must always be discarded after use. Re-use of an irrigation line or clip may result in contamination of your patients. Do not attempt to resterilize the irrigation lines. Do not attempt to modify the irrigation lines and clips.

10.2 CONTRA-ANGLE MAINTENANCE

Refer to the manufacturer's instructions for your contra-angle.

10.3 ROTARY INSTRUMENT MAINTENANCE

Please refer to the manufacturer's instructions for your instruments (burs, drills, etc.).

10.4 DEVICE MAINTENANCE

Important:

Do not use an abrasive product to clean the device.

Do not allow cleaning/disinfecting spray to penetrate inside the device.

The device, the bracket and the multifunction footswitch are not sterilizable.

Unplug the main unit from power source before cleaning or disinfecting. Do not use directly on liquids and sprays. Never immerse the unit or its parts in liquid.

The housing is not waterproof. Do not use ultrasonic or steam sterilizers.

The HT IMPLANT control unit, the control footswitch and the bracket must always be cleaned and disinfected after each procedure using alcohol, disinfectant agents or disinfectant wipes for dental practices.

Regular monitoring of the HT IMPLANT unit is necessary in order to detect any problem.

10.5 MICRO-MOTOR AND CORD MAINTENANCE

Do not sterilize or lubricate the micro motor. Do not attempt to disassemble the motor or motor connector. The saline solution corrodes all the metal parts. Immediately remove residue.

External surface of the micro motor may be cleaned with a damp cloth or disinfectant wipe. Do not place micro motor in ultrasonic sterilizers or submerge in any solutions.

10.6 MICRO-MOTOR HOLDER MAINTENANCE

Cleaning and disinfection

Motor Holder may be cleaned with a damp cloth or disinfectant wipe.

When treating patients who may have an acute, critical infectious disease, be sure to observe the hygienic measures cited in applicable publications and reports. If possible, use suitable disposable products to avoid the transmission of critical pathogens. These protect the user, the patient and all participants in the surgery. A microbiologically appropriate coolant liquid must be used. Only use prescribed isotonic saline solution NaCl 0.9 (identified as an infusion solution; follow the instructions on the package insert) for cooling and spraying wounds.

XI - MONITORING / PREVENTIVE AND CORRECTIVE MAINTENANCE

Important:

In the case of a fault, you are advised to contact the supplier of your device. Do not use unauthorized repairers, who might make your device dangerous for you and your patients.

11.1 MONITORING

Regular monitoring of the device and its accessories is necessary in order to detect any insulation fault or damage. Replace, if necessary.

It is important to keep the control unit ventilation vents clean in order to avoid abnormal heating.

11.2 PREVENTIVE AND CORRECTIVE MAINTENANCE

HT IMPLANT micro-motor inspection/overhaul: The manufacturer recommends inspection or overhaul of the HT IMPLANT micro-motor at least once a year.

The micro motor cable has limited durability and is considered a consumable product. The durability of the cable is largely dependent upon use and treatment, i.e., handling, reprocessing and frequency of use. Caution must be exercised when connecting or disconnecting the micro motor and foot control cable from the device to avoid damage to the cable (Fig. 8).

Periodically inspect the micro motor cord for damage prior to use. If there is damage noted, contact SATELEC.

1. Do not attempt to disassemble the motor or motor connector.
2. Do not oil or lubricate the motor.
3. Do not attach or remove a handpiece to the motor while the motor is running.
4. Do not bend motor cord sharply.
5. Do not drop or impact micro motor against hard surface.

Failure to comply with any of the above instructions may void your warranty

The average life of this medical device with normal use and in accord with the present user manual, has been calculated to be 5 years. According with EN62353 we recommend users verify every 2 years from the date of the first use: dispersion current versus ground, in the box and in the patient.

MOTOR RUBBER O-RING:

- If necessary, substitute the "Rubber O-Ring"

LUBRICATING:

- **Do not oil or lubricate the micro motor.** Protect motor from excess oil draining from handpiece. After lubricating and before autoclaving, stand handpiece by its base on a paper towel and allow excess oil to drain.
- The internal gears of the peristaltic pump are self lubricating.
- Do apply some oil drops in the cylinders of the peristaltic pump rotor after a long period of non-use or when the pump doesn't run continuously.

11.3 FUSE REPLACEMENT

Your device is protected by two fuses located in the power receptacle. The replacement procedure is as follows (Fig. 11):

- **ONCE UNIT IS UNPLUGGED** -With a screwdriver remove the fuse holder from the power inlet connector.
- Extract the fuse holder, replace the fuse with one of the same type and value 2 X 1A T(230V) / 2 X 2A T(115V) (Dim. 5x20).
- After inspection or replacement of a fuse, re-secure fuse holder ensure that the correct voltage value appears.

Note:

All the information necessary for the repair of faulty components is made available by the manufacturer to the technical personnel of SATELEC approved dealer network who carry out the work.

- The HT IMPLANT must be used in accordance with these Instructions for Use.
- The HT IMPLANT has no components which can be repaired by the user. Assembly, modifications or repairs must only be undertaken by an authorized service organization.
- Unauthorized opening of the equipment invalidates all claims under warranty and any other claims.
- Use on electrical circuit complying with the terms of regulation CEI 64-8 section 710.

11.4 OPERATION FAULTS

Fault detected: The central unit does not switch on.

Possible cause	Solution
• Mains cord not connected properly	Make sure the Mains cord is properly plugged into the electrical outlet and the connector on the back of the unit.
• Unit has not been switched on.	Check the back of the unit to see that the On button has been pressed.
• Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another working electrical outlet. If the unit does not turn on, verify the integrity of the fuses in the unit.
• Broken connections inside. • Transformer damaged. • PCB faulty. • Filter supply defective. • Power interrupted	Immediately discontinue use of the device and contact SATELEC after-sales service.

Fault detected: The central unit turns off and stops working suddenly

Possible cause	Solution
• Mains cord damaged. • Filter damaged. • Transformer damaged. • PCB damaged. • Connections inside interrupted.	After checking the correct connection to the electrical outlet and the absence of liquid in the unit, try to power the central unit. If the device does not turn on, immediately discontinue use of the device and contact SATELEC after-sales service.

Fault detected: The central unit doesn't switch off

Possible cause	Solution
• Power filter damaged.	Unplug the mains cord from the electrical outlet and the central unit. Then contact SATELEC after-sales service.

Fault detected: The footswitch doesn't start the motor

Possible cause	Solution
• Micro Motor cable not connected properly	Check to see that the Micro Motor cable is securely plugged into the connector on the front right of the central unit, disconnect and reconnect.
• Footswitch cable is not connected properly	Check that the footswitch cord is plugged into the connector on the front left of the central unit, disconnect and reconnect.
• Male connector of the Micro motor damaged.	Unplug the Micro Motor and make sure that the male connector is not damaged, wet or damp.
• Pins of the connector damaged	Unplug the Micro Motor and the footswitch cable, make sure the metal pins on the male connector are intact and not damaged and have not broken off and/or remained stuck in the female connector mounted on the unit.
• Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.
• Footswitch damaged.	Check the ability to start the Micro motor from the central unit keyboard, pressing the [Motor ON/OFF switch]. If so, contact SATELEC after-sales service.

Fault detected: It is not possible to change “Program”, FWD/REV and Flow by using the footswitch.

Possible cause	Solution
• Footswitch damaged. • Metal pins damaged in the connector	Check the possibility of changing the parameters by using the unit keyboard. Contact SATELEC after-sales service.

Fault detected: The display is not working properly

Possible cause	Solution
• Contra-angle damaged.	Check the integrity of the Contra-angle.
• Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.
• Possible interference.	Check to see if the central unit has been positioned too close to another electrical device which can cause interference. In this case, move the interfering electrical device further away from the central unit.
• Broken connections inside. • Display damaged. • PCB damaged.	If you cannot find the cause of the malfunction, immediately discontinue use of the device and contact SATELEC after-sales service.

Fault detected: The keyboard is not working properly

Possible cause	Solution
• Keyboard damaged. • Broken connections inside.	Attempt to change parameters by using the footswitch. Contact SATELEC after-sales service.

Fault detected: The sound signals are not working properly

Possible cause	Solution
• Broken connections inside	Immediately discontinue use of the device and contact SATELEC after-sales service.

Fault detected: The motor doesn't work

Possible cause	Solution
• Micro motor cable not connected properly	Check to see that the micro motor cable is plugged into the connector on the front right of the central unit, disconnect and reconnect.
• Contra-angle defective or not properly calibrated.	Check the integrity of the Contra-angle. Remove the Contra-angle and perform calibration.
• Footswitch damaged.	Attempt to start the Micro Motor by using the [Motor ON/OFF button] on the central unit keyboard. If unit does not turn on, contact SATELEC after-sales service.
• Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.
• PCB damaged. • Motor damaged. • Connector pins and/or cable damaged.	Immediately discontinue use of the device and contact SATELEC after-sales service.

Fault detected: The motor loses power.	
Possible cause	Solution
• Contra-angle defective or not properly calibrated.	Check the integrity of the Contra-angle by removing the Contra-angle and perform calibration.
• Motor damaged.	Contact SATELEC after-sales service.
Fault detected: The motor stops suddenly	
Possible cause	Solution
• Micro motor/Cable damaged	Check to see that the micro motor cable is plugged into the connector on the front right of the central unit, disconnect and reconnect. Start the motor by pressing the [Motor ON/OFF] button on the keyboard. If the motor turn on, immediately discontinue use of the device and contact SATELEC after-sales service.
• Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.
• Footswitch damaged.	Attempt to start the Micro Motor by using the [Motor ON/OFF button] on the central unit keyboard. If unit does not turn on, contact SATELEC after-sales service.
• PCB damaged.	If you cannot find the cause of the malfunction, immediately discontinue use of the device and contact SATELEC after-sales service.
Fault detected: The Motor will not turn off	
Possible cause	Solution
• Keyboard and/or footswitch damaged	Try to turn off the motor by pressing the [Motor ON/OFF] button on the keyboard or by using the foot pedal. If the motor will not turn off, immediately discontinue use of the device and contact SATELEC after-sales service.
• Broken connections inside. • PCB damaged	Immediately discontinue use of the device and contact SATELEC after-sales service.
• Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.
Fault detected: The pump is not working property	
Possible cause	Solution
• Rotor is blocked.	Disconnect power to the central unit. Check to see if the rotor of the peristaltic pump rotates correctly and is not blocked by any internal or external objects. Remove any visible blockage if possible. Otherwise, contact SATELEC after-sales service.
• Tube not correctly inserted • Tube damaged.	To verify the proper tube connection follow the instructions on Chapter VIII - USING THE DEVICE. Check for damage along the tube, kinks, partial occlusions, etc. Arrange for immediate replacement of the tube.

Fault detected: The pump stops suddenly

Possible cause	Solution
•Rotor is blocked.	Disconnect power to the central unit. Check to see if the rotor of the peristaltic pump rotates correctly and is not blocked by any internal or external objects. Remove any visible blockage if possible. Otherwise, contact SATELEC after-sales service.
•Pump damaged/blocked. •PCB damaged.	Immediately discontinue use of the device and contact SATELEC after-sales service.

XII - ELECTROMAGNETIC COMPATIBILITY

Important:

The power cord, micro-motor cord, and footswitch cord must be kept apart.

The HT IMPLANT requires special precautions to be taken with regard to electromagnetic compatibility. It must be installed and prepared for use as described in chapter IV.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the HT IMPLANT. The recommended separation distances in this paragraph must therefore be complied with.

The HT IMPLANT must not be used near to or on top of another device. If this cannot be avoided, its operation under the conditions of use must be checked beforehand.

The use of accessories other than those specified or sold by BY DENTAL and/or distributed by SATELEC as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the HT IMPLANT.

To avoid possible risks of electromagnetic interferences do not use other electro-medical instruments of similar nature near HT IMPLANT. The unit complies with the current normative electromagnetic radiation law.

This device has been tested and found to comply with the emissions requirements of IEC 60601-1-2. These requirements provide reasonable protection against harmful electromagnetic interference in a typical medical installation. However, high levels of radio-frequency (RF) emissions from electrical devices, such as cellular phones, may disrupt the performance of this device. To mitigate disruptive electromagnetic interference, position this device away from RF transmitters and other sources of electromagnetic energy.

Risks from electromagnetic fields: The functions of implanted systems (such as pacemakers) can be influenced by electromagnetic fields. Patients and users must inquire about implanted systems before using the motors and check the use. Do not bring the motors close to the systems. Take suitable emergency precautions and immediately react to any changes in health.

12.1 ELECTROMAGNETIC EMISSIONS

HT IMPLANT is intended for use in the electromagnetic environment specified below. The HT IMPLANT should only be used in such an environment.

RF emissions CISPR 11	Group 1	HT IMPLANT uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	HT IMPLANT is suitable for use in all establishment including domestic establishments and those directly connected to the public low-voltage power supply.
Harmonic emissions IEC 61000-3-2	Class A Conforms	
Voltage fluctuation/ flicker emissions IEC 61000-3-3	Conforms	

12.2 ELECTROMAGNETIC IMMUNITY

HT IMPLANT is intended for use in the electromagnetic environment specified below. The HT IMPLANT should only be used in such an environment.

Immunity test	Test level IEC 60601-1-2	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) EN61000-4-2	±6kV contact ±8kV air	±6kV ±8kV	Floor should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%
Burst/Fast transient EN 61000-4-4	±2kV for power supply lines	±2kV for power supply lines	Quality of the Main power source should be that of a typical commercial or hospital environment
Surge EN 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Quality of the Main power source should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variation on power supply input lines. EN 61100-4-11	<5%U _t (>95% dip in U _t) for 0,5 cycle 40% U _t (60% dip in U _t) for 5 cycle 70% U _t (30% dip in U _t) per 25 cycle	<5%U _t (>95% dip in U _t) for 0,5 cycle 40% U _t (60% dip in U _t) for 5 cycle 70% U _t (30% dip in U _t) per 25 cycle	Quality of the Main power source should be that of a typical commercial or hospital environment If the user of the HT IMPLANT requires continued operation during power main interruptions, it is recommended than the HT IMPLANT be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field EN 61100-4-8	3A/m	3A/m	Power frequency magnetic fields should be at a typical commercial or hospital establishment level.

HT IMPLANT is intended for use in the electromagnetic environment specified below. The HT IMPLANT should only be used in such an environment.

Immunity test	Test level IEC 60601-1-2	Compliance level	Electromagnetic Environment - Guidance
Conducted RF EN 61000-4-6	3 Veff from 150KHz to 80MHz	3 Veff from 150KHz to 80MHz	Portable and mobile RF communication equipment should be used no closer to any part of the HT IMPLANT, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d = 1,2√P from 150KHz to 80MHz d = 1,2√P from 80MHz to 800MHz d = 2,3√P from 800MHz to 2.5GHz where P is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer and recommended separation distance in metres (m).
Radiated RF IEC 61000-4-3	3 Veff from 80MHz to 2.5GHz	3 Veff from 80MHz to 2.5GHz	
Field strengths from fixed RF transmitter, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b . Interference may occur in the proximity of equipment marked with the following symbol: 			

12.3 RECOMMENDED SEPARATION DISTANCES

Recommended separation distance between portable and mobile RF communications equipment and the HT IMPLANT			
The HT IMPLANT is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the HT IMPLANT can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the HT IMPLANT as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according in frequency of transmitter (m)		
	From 150kHz to 80MHz $d = 1,2\sqrt{P}$	From 80MHz to 800MHz $d = 1,2\sqrt{P}$	From 800MHz to 2,5GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance (D) in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

XIII - DISPOSAL AND RECYCLING

As electrical and electronic equipment, the device must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (in particular on the European market, with reference to Directive 2002/96/EC of 23/01/2003).

When your device reaches the end of its life, we consequently recommend that you contact your dental equipment dealer (or, failing this, the nearest ACTEON GROUP office, the list of which is given in chapter XVIII), for information on how to proceed.

XIV - MANUFACTURER'S LIABILITY

The manufacturer is not liable if:

- the manufacturer's installation recommendations have not been followed (supply voltage, electromagnetic environment, etc.),
- repairs or other work have been done on the device by persons not authorized by the manufacturer,
- the device has been used in an electrical installation which does not comply with current standards,
- the device has been used in a way which is not stipulated in this manual,
- accessories (motor and cord, irrigation line, etc.) other than those supplied by BY DENTAL and/or distributed by SATELEC have been used,
- the instructions in this document have not been followed.

Note:

The manufacturer reserves the right to modify the device and/or the user manual without notice.

XV - REGULATIONS

This medical device is classified as class IIa according to European directive 93/42/EEC. This equipment is manufactured in compliance with the current IEC 60601-1 standard. This equipment has been designed and manufactured according to an ISO 13485-certified quality assurance system.

XVI - SYMBOLS

This **Operating and Maintenance Instruction Manual** provides instructions on HT IMPLANT. The text, graphics and icons are used to facilitate understanding of the functions and operations by the final user. Target Group: Dentists and office personnel.

It is very important to keep this manual with the device, for easy reference at any time in the future. In the event of sale or transfer to another final user, be sure to attach this Operating and Maintenance Instruction Manual. Transferring this manual will ensure, the new final user will have all the necessary information to understand the functionality and proper use of the unit. These instructions are intended to ensure the safe use of this equipment and should be read carefully before installation and use of the device.

Description of safety instruction


	Important information, explanation or integration for users.
	Have a direct impact on the functionality of this product or warns of possible misuse and/or harm to the product or the user.

Three hazard levels are used in this document for avoiding personal and property damage:

	Indicates a hazardous situation that can cause damage to property, or mild or moderate physical harm. (*)
	Indicates a hazardous situation that can cause death or serious injury. (*)
	Indicates the maximum hazard level. Indicates a directly hazardous situation that can cause death or serious injury. (*)

(*) The optional step covers necessary measures for avoiding hazards.

Symbols

	EEC 93/42		Risk of electric shock
	Follow Instructions for use		Disposal Information
	Type BF appliance	V	Supply Voltage
	Manufacturer	REF	Reference
		Type	Device type
SN	YY = year manufactured ZZZ = Ref XXXX = Serial Number		Installation side and securing the mains input
	Fuse		Mode: continuous operation with intermittent load
	No access for persons with active implant as pacemaker		

The symbols printed on the outside are for transportation and storage, and have the following meaning:

	Protect from moisture		Transport upright with the arrows pointing upwards
	Fragile - Protect against impact!	Handle With Care	=

	<i>Please keep the packaging in case you need to return the product for servicing or repair.</i>
---	--

Refer to the [SN=Serial Number](#) in any communication with the Manufacturer.

This symbol placed on the above sticker remind to consult the present Operating and Maintenance Instruction Manual

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